

# Investigating Injuries: Patellar Tendinopathy



# Introductions



Dr. Daniel G. Stewart  
PT, DPT

 [dr.danielstewartdpt](https://www.instagram.com/dr.danielstewartdpt)






Michael Morgan  
PT, DPT, OCS, CP-MSK

 michael-morgan-dpt-ocs

 Michael-Morgan-60



# Learning Objectives




- Describe the anatomy, biomechanics, and pathophysiology of Patellar Tendinopathy
- Perform a structured clinical assessment for Patellar Tendinopathy
- Utilize ActivForce handheld dynamometry to objectively assess quadriceps strength and tendon loading capacity
- Explain the clinical role of musculoskeletal ultrasound imaging and Blood Flow Restriction (BFR) training

# Agenda



- Introduction: Understanding Patellar Tendinopathy
- Clinical Presentations
- Mechanisms of injury
- Common Causes and Risk Factors
- Clinical Assessment
  - Range of Motion (ROM)
    - Knee Flexion/Extension
  - Muscle Testing
    - Quadriceps Strength
    - Hip Abductors and Extensors
- Case Study
  - Smartphone Clinicians Guide to Patellar Tendinopathy: Use of ActivForce, Ultrasound Imaging and Blood Flow Restriction
  - Data Review and Analysis
- Treatment and Rehabilitation Strategies
- Q&A

# Understanding Patellar Tendinopathy





# What is Patellar Tendinopathy?

- Also known as “Jumper’s Knee”
- Overuse injury affecting the knee extensor mechanism
- Common in jumping and explosive sports
- Characterized by localized pain at the inferior pole of the patella
- Typically associated with load-related pain





# What is Patellar Tendinopathy?

- Common in athletes aged 15–30
- Higher prevalence in elite jumping athletes
- Reported prevalence:
  - ~45% in elite volleyball athletes
  - ~14% in recreational athletes
- More common in male athletes



# Clinical Presentations



## Common symptoms:

- Localized pain at inferior pole of patella
- Pain with jumping, squatting, or running
- Morning stiffness
- Pain during resisted knee extension
- Reduced athletic performance

# Anatomy



Structures involved in the knee extensor mechanism:

- Quadriceps muscle group
- Quadriceps tendon
- Patella (sesamoid bone)
- Patellar tendon
- Tibial tuberosity

Function:

- Transmission of force from quadriceps to tibia
- Enables knee extension
- Important for shock absorption during landing






# Anatomy - Quadriceps Muscle Group

## Primary knee extensor muscles

- Rectus femoris
- Vastus lateralis
- Vastus medialis
- Vastus intermedius

## Function

- Knee extension
- Eccentric control during landing
- Shock absorption during deceleration



# Anatomy - Quadriceps Muscle Group



Link: [youtube.com/shorts/kL1Y-D4P5Q4](https://www.youtube.com/shorts/kL1Y-D4P5Q4)




# Anatomy - Patellar Tendon Structure

- Connects inferior patella to tibial tuberosity
- Composed primarily of Type I collagen fibers
- Designed to handle high tensile loads
- Works as part of the extensor mechanism

## Mechanical role

- Transfers quadriceps force to extend the knee
- Stores and releases elastic energy during jumping





# Mechanism of Injury



- Chronic mechanical overload
- Repetitive micro-tearing of tendon fibers
- Degenerative changes:
  - Collagen disorganization
  - Increased ground substance
  - Neovascularization

## Common mechanisms:

- Repetitive jumping and landing
- Sudden increase in training volume
- Repeated knee extension loading
- Insufficient recovery between sessions



# Mechanism of Injury

## Sports commonly associated:

- Volleyball
- Basketball
- Track and field
- Soccer



Link: [youtube.com/shorts/7\\_ErV4C-hcs](https://www.youtube.com/shorts/7_ErV4C-hcs)

**Key Concept:** Tendinopathy is degenerative (tendinosis) rather than inflammatory.



# Risk Factors

## Intrinsic Factors

- Quadriceps tightness
- Hamstring tightness
- Increased Q-angle
- Muscle imbalance
- Reduced ankle dorsiflexion

## Extrinsic Factors

- Sudden increase in training
- High training volume
- Hard playing surfaces
- Poor biomechanics

# Clinical Assessment:



# Range of Motion



## Knee

- Flexion
- Extension



# Range of Motion



## Hip

- Extension
- Rotation



# Muscle Testing



## Quadriceps

- Knee Extension Strength



# Muscle Testing

## Hip Muscles

- Hip abductors
- Hip extensors



# Case Study

A Data-Driven Approach to Anterior Knee Pain:  
Integrating Dynamometry, Ultrasound, and BFR





# A tale of two presentations of inferior patellar pain

## Case 1

- Male
- Mid 30s
- Distance trail runner
  - Increased frequency of running Camelback Mountain in Phoenix, AZ
  - Training for Rim2Rim2Rim Grand Canyon run (~45 miles/10k elevation)
- 3 month history of worsening infrapatellar pain



## Case 2

- Female
- Mid 30s
- Distance trail runner
  - Increased frequency & duration of indoor cycling in summer due to outdoor heat
- 1 month history of worsening infrapatellar pain

# Objective Exam



Exam Finding	Case 1	Case 2
Primary pain location (patient report)	Proximal patellar tendon (PT)	Proximal patellar tendon (PT)
Active knee flexion	Full and pain-free	Sharp end range pain
Active knee extension	Full and pain-free	Sharp end range pain
Selective Tissue Testing: Knee Extensors	<b>Strong and painful</b>	<b>Strong and pain free</b>
Palpation	TTP over proximal PT - mild with deep pressure	TTP over proximal PT- exquisite with deep pressure
Ligamentous/Meniscal Testing	Negative	Negative
Step down test	Reproduces anterior knee pain localized to tendon	Minimal symptoms



# Ceiling effect of MMT for Knee Extension

> [Percept Mot Skills](#). 2000 Jun;90(3 Pt 1):747-50. doi: 10.2466/pms.2000.90.3.747.

## A broad range of forces is encompassed by the maximum manual muscle test grade of five

R W Bohannon<sup>1</sup>, D Corrigan

Affiliations + expand

PMID: 10883753 DOI: [10.2466/pms.2000.90.3.747](https://doi.org/10.2466/pms.2000.90.3.747) [↗](#)

### Abstract

This retrospective study summarizes the range of knee extension forces associated with the manual muscle test grade of 5. The range of the forces (>560 Newtons) was large and encompassing of more than 86% of measurable forces. These findings help to explain the insensitivity of manual muscle testing and why it has such a profound ceiling effect.


# Dynamometry



## Case 1 Initial Evaluation

### Knee Extension Seated

#### Force over time



#### Average Values

Right	27.58 lb
Left	44.31 lb


#### Peak Force (lb)

Right	36.02 lb
Left	53.11 lb
Strength Difference	17.08 lb
Percentage Difference	38.35%

## Case 2 Initial Evaluation

### Knee Extension Seated

#### Force over time




#### Average Values

Right	35.13 lb
Left	34.70 lb

#### Peak Force (lb)



Right	39.68 lb
Left	39.33 lb
Strength Difference	0.35 lb
Percentage Difference	0.89%

# Ultrasound Imaging





Tendon State	Thickness
Mild tendinopathy	4–5 mm
Moderate tendinopathy	5–7 mm
Severe tendinopathy	>7 mm

# Case 1 - Ultrasound Confirmed Patellar Tendinopathy



**Focal proximal tendon thickening (~5-6 mm compared to 4mm contralaterally) with normal Doppler activity supports patellar tendinopathy rather than fat pad irritation.**

# Case 2 - Ultrasound Confirmed Fat Pad Syndrome



Rheumatol Int (2017) 37:757–764  
DOI 10.1007/s00296-016-3647-4

Rheumatology  CrossMark  
INTERNATIONAL

IMAGING

**Sonographic characterization of Hoffa's fat pad. A pilot study**

Erika Vera-Pérez<sup>1</sup> · Guadalupe Sánchez-Bringas<sup>2,3</sup> · Lucio Ventura-Ríos<sup>3</sup> ·  
Cristina Hernández-Díaz<sup>3</sup> · Socorro Cortés<sup>4</sup> · Marwin Gutiérrez<sup>3</sup> · Carlos Pineda<sup>5</sup>

**Findings consistent with hypervascularity of the infrapatellar fat pad without evidence of patellar tendinopathy.**

# Treatment & Rehabilitation Strategies



# Treatment



	<b>Case 1 – Patellar Tendinopathy</b>	<b>Case 2 – Hoffa’s Fat Pad Irritation</b>
Primary pathology	Proximal patellar tendon overload	Infrapatellar fat pad irritation
Strength testing	<b>Extension strength deficit on dynamometry</b>	<b>Symmetrical strength bilaterally</b>
Initial symptom driver	Load intolerance	Mechanical compression
<b>Primary management focus</b>	Progressive tendon loading	Mechanical modification
Activity modification	Reduce high-load activities (running, jumping)	Temporary cycling modification
Key intervention	<b>BFR-assisted strengthening protocol</b>	<b>Bike fit correction + anti-inflammatory care</b>



# Case 2 - Hoffa's Fat Pad Management

Intervention	Strategy
Mechanical correction	Crank arm length reduced <b>175 mm</b> → <b>160 mm</b>
Saddle position	Adjusted to <b>heel-on-pedal method</b> for proper knee extension
Load modification	Temporary reduction in cycling intensity
Anti-inflammatory treatment	<b>Iontophoresis with dexamethasone</b> prescribed by PCP
Outcome	<b>Symptoms resolved within one week with no recurrence</b>




Pain was driven by **mechanical compression rather than contractile deficit**, so correction of biomechanics resolved symptoms rapidly.



# Case 1 - Intervention



Phase	Intervention
<b>Weeks 0–4</b>	Blood Flow Restriction strengthening <b>2–3x/week</b>
	Load prescribed using <b>dynamometry-derived % max force</b>
	Low-load quad focused strengthening (20–40% 1RM) @ 70% LOP
<b>Weeks 4–8</b>	Transition to <b>progressive heavy slow resistance program at gym</b>
<b>Weeks 8-12</b>	<b>Power/plyometrics</b> incorporated into gym-based strengthening
<b>Outcome</b>	<b>Returned to hiking pain-free at 12 weeks &amp; running at 16 weeks</b>




Isometric Knee Extension + BFR	Open Chain Knee Extension + BFR	Closed Chain Step-Down + BFR
<p>Establish <b>tolerance + force output</b></p> <ul style="list-style-type: none"> <li>Use <b>ActivForce</b> to determine %MVC &amp; provide <b>biofeedback</b></li> </ul>	<ul style="list-style-type: none"> <li>Load based on prior <b>MVC output</b></li> <li>Progress to <b>dynamic contractions</b></li> </ul>	<ul style="list-style-type: none"> <li>Transition to <b>functional loading</b></li> <li>Emphasize <b>eccentric control + alignment</b></li> </ul>

# Case 1 - Outcome

## Knee Extension Seated



Force over time



### Average Values

Right	50.35 lb
Left	48.48 lb

### Peak Force (lb)

Right	58.01 lb
Left	55.56 lb
Strength Difference	2.46 lb
Percentage Difference	4.31%





# Data Driven Approach to Patellar Tendon Pain

- **Measure the Impairment**
  - Dynamometry provides objective strength assessment beyond manual muscle testing
- **Identify the Tissue**
  - Diagnostic ultrasound helps differentiate tendon pathology from adjacent structures such as Hoffa's fat pad
- **Target the Treatment**
  - Tendinopathy → progressive loading strategies (BFR when load intolerance exists)
  - Fat pad irritation → biomechanical correction and symptom management
- **Objective data improves clinical decision making.**

Dynamometry and ultrasound together help physical therapists refine diagnosis and guide more precise treatment strategies.
- **Better measurement → better diagnosis → better rehabilitation.**

# Questions?

SMART  CUFFS

 **clarius** Ultrasound  
Anywhere



**POCUS**

Point-of-Care Ultrasound  
Certification Academy™





# References

- Cook, J. L., & Purdam, C. R. (2009). Is tendon pathology a continuum? *British Journal of Sports Medicine*.
- Malliaras, P., Cook, J., & Purdam, C. (2015). Patellar tendinopathy: Clinical diagnosis and management. *British Journal of Sports Medicine*.
- Walton, D., et al. (2023). Patellar tendinopathy in athletes involved in jumping sports. *Sports Health*.
- StatPearls Publishing. (2024). Patellar Tendinopathy. *NCBI Bookshelf*.
- Mayo Clinic. (2024). Patellar Tendinitis (Jumper's Knee).
- TeachMe Orthopedics. (2024). Jumper's Knee (Patellar Tendinopathy).
- Vaughan Physiotherapy. (2024). Patellar Tendinopathy / Jumper's Knee.
- Musculoskeletal Key. (2023). Pathophysiology of Tendinopathy.